## **Platelet Function Studies**

Order Name: PLT FUN
Test Number: 1506325
Revision Date: 08/26/2014

TEST NAME	METHODOLOGY	LOINC CODE
Platelet Function, ADP	Platelet Function Testing	24472-3
Platelet Function, Epinephrine	Platelet Function Testing	24471-5

SPECIMEN REQUIREMENTS						
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment		
Preferred	5.4 mL	Whole Blood	Sodium Citrate 3.2% (Blue Top)	Room Temperature		
Instructions	NOTE: If collected at a location other than the laboratory at St. John Medical Center 1923 South Utica Ave. Tulsa, Then please send by STAT courier the the laboratory at St. John Medical Center for testing.  Specimen must be tested within 4 hours of collection. Do not refrigerate!  Collect Two 2.7mL Sodium Citrate 3.2% (Blue Top) tubes - DO NOT Spin, Filter or Freeze specimens!					
	Patient should have PLT >150  Please write on request if pa	,000 and HCT >35% for accuraction is receiving aspirin.	cy.			

GENERAL INFORMATION		
Testing Schedule	Daily	
Expected TAT	1 Day	
Clinical Use	Platelet function studies are done to evaluate platelet function. This is a specialized test and would normally be performed in patients with some indicator of a qualitative platelet disorder.	
CPT Code(s)	85576x2	
Lab Section	Coagulation	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025