Labcorp Oklahoma, Inc. Test Directory

Ristocetin Platelet Aggregation (High and Low)

Order Name: **RIPA AGG** Test Number: 1501930 Revision Date: 10/25/2023

METHODOLOGY	
Aggregation	5998-0
Aggregation	
Flow cytometry	777-3
	Aggregation Aggregation

Pathology Report

SPECIMEN REQUIREMENTS						
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment		
Preferred	See Instructions	Whole Blood	Sodium Citrate 3.2% (Blue Top) and EDTA (Lavender Top)	Room Temperature		
Instructions	READ BEFORE COLLECTING SPECIMEN!					
	COLLECTIONS and TESTING MUST BE SCHEDULED with the lab at Ascension St. John Coagulation Department! Please call the coagulation department to make testing arrangements: (918) 744-2500 x15513. If OFF-SITE collection is Authorized by the Coagulations Department, the specimens must reach Ascension St. John Main Lab within 1 hour of collection and before 1pm. If you cannot arrange for specimens to arrive in this time frame, Do Not Collect Specimen! Patient Must be fasting for at least 8 hours before collection.					
	Collect Both: - Four (2.7mL) 3.2% Sodium Citrate - Blue top tubes - One (4.5mL) EDTA - Lavender top tubes					
	Keep specimens whole blood (Do Not Spin) - Keep specimen at Room Temperature!					
	Patient should refrain from aspirin, Greater than 75,000 for accuracy.	phenylbutazone, phenothiazines or	antihistamines for 10 days prior to the test. F	Patient should have PLT count		

GENERAL INFORMATION		
Testing Schedule	Mon - Fri (7am - 1pm)	
Expected TAT	1 Day	
Clinical Use	Platelet aggregation studies are done to evaluate platelet function. This is a specialized test and would normally be performed in patients with some indicator of a qualitative platelet disorder.	
CPT Code(s)	85576x2, 8557626	
Lab Section	Coagulation	