## Labcorp Oklahoma, Inc. Test Directory

METHODOL OGV

## Serotonin, Whole Blood

Order Name: SEROTON B
Test Number: 3602470
Revision Date: 12/12/2022

| IESI NAME            |   |                  | METHODOLOGY                          | LOING CODE            |
|----------------------|---|------------------|--------------------------------------|-----------------------|
| Serotonin, Whole Blo | od  |                  | Liquid Chromatography/Tandem Mass Sp | ectrometry 2939-7     |
| SPECIMEN REQUIRE     | EMENTS  |                  |                                      |                       |
| Specimen             | Specimen Volume (min)   | Specimen Type    | Specimen Container                   | Transport Environment |
| Preferred            | See Instructions  | See Instructions | LabCorp ascorbic acid transpo        | ort Frozen            |
| Instructions         | Specimen Type: Lavender-top (EDTA) tube and transport tube (LabCorp item ID 33314) containing 75 mg ascorbic acid Specimen Storage: FREEZE Please note that stability claims reflect stability of properly prepared whole blood with ascorbic acid preservative.  Specimen Collection: Sample collection practices that prevent blood clotting and oxidation of serotonin are imperative. Collect three tubes (4 mL each) whole blood. Invert tubes at least ten times to ensure thorough distribution of EDTA. IMMEDIATELY transfer 10 mL whole blood to transport tube containing 75 mg ascorbic acid. Thoroughly mix by inverting whole blood/ascorbic acid mixture another ten times. FREEZE transport tube and SEND ENTIRE COLLECTION to the laboratory. To avoid delays in turnaround time when requesting multiple tests on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.  Specimen Stability: Ambient: 4 hours, Refrigerated: 24 hours, Frozen: 16 days  Medications that may affect serotonin concentrations include lithium, MAO inhibitors, methyldopa, morphine, and reserpine. In general, foods that contain serotonin do not interfere significantly. Slight increases may be seen in acute intestinal obstruction, acute MI, cystic fibrosis, dumping syndromes, and nontropical sprue. Metastasizing abdominal carcinoid tumors often show serotonin concentrations greater than 400 ng/mL. |                  |                                      |                       |

| GENERAL INFORMATION             |  |  |
|---------------------------------|--|--|
| Expected TAT                    | 3 - 6 days   |  |
| Clinical Use                    | This is the preferred serotonin test when diagnosing carcinoid tumors. |  |
| Performing Labcorp Test<br>Code | 120089   |  |
| Notes                           | Labcorp Test Code: 120089  |  |
| CPT Code(s)                     | 84260  |  |
| Lab Section                     | Reference Lab  |  |
|                                 |  |  |

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025