## Labcorp Oklahoma, Inc. Test Directory

## Inhibin A

Order Name: INHIBIN A
Test Number: 3622375
Revision Date: 12/12/2022

TEST NAME	nE		DDOLOGY	LOINC CODE	
Inhibin A		Immun	ochemiluminometric (ICMA)	23883-2	
SPECIMEN REQUIREM	MENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	0.3 mL (0.2)	Serum	Clot Activator SST	Frozen	
Alternate 1	0.3 mL (0.2)	Serum	Clot Activator (Red Top, No-Gel)	Frozen	
Instructions	Notes: 0.2 mL (Note: This volume Does NOT allow for repeat testing.)  Specimen Type: Red-top tube or gel-barrier tube  Specimen Storage: Freeze  Specimen Collection: Transfer the serum into a LabCorp PP transpak frozen purple tube with screw cap (LabCorp ID 23751). The specimen should be FROZEN immediately and maintained frozen until tested. To avoid delays in turnaround time when requesting multiple tests on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.				
	<b>Special Instructions:</b> This test may exhibit interference when sample is collected from a person who is consuming a supplement with a high dose biotin (also termed as vitamin B7 or B8, vitamin H, or coenzyme R). It is recommended to ask all patients who may be indicated for				
	,	pecimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available			
	Specimen Stability. Ambient, Not Available, Nethigerated : Not Available, 1 102611. Not Available				

GENERAL INFORMATION		
Expected TAT	4 - 6 days	
Performing Labcorp Test Code	146803	
Notes	Labcorp Test Code: 146803	
CPT Code(s)	86336	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025