## Heavy Metals Panel, Urine (24hr or Random) with Reflex to Arsenic Fractions

Order Name: METAL S U
Test Number: 3810575
Revision Date: 11/15/2023

TEST NAME METHODOLOGY LOINC CODE

Heavy Metals Panel, Urine (24hr or Random) with Reflex to Arsenic Fractions

Atomic Absorption Spectrophotometry

| SPECIMEN REQUIREMENTS |   |                |   |                       |  |  |
|-----------------------|---|----------------|---|-----------------------|--|--|
| Specimen              | Specimen Volume (min)   | Specimen Type  | Specimen Container                            | Transport Environment |  |  |
| Preferred             | 15 mL (5.3 mL)  | Urine, 24-hour | 24 Hour Urine Acid Washed<br>Container        | Room Temperature      |  |  |
| Alternate 1           | 15 mL (5.3 mL)  | Urine, Random  | Acid Washed, Trace Element Free<br>Contatiner | Room Temperature      |  |  |
| Instructions          | Notes: 5.3 mL (Note: This volume Does NOT allow for repeat testing.)  Specimen Type: Plastic urine container, no preservative  Specimen Storage: Maintain specimen at room temperature.  Specimen Collection: Sampling time is not critical for industrial monitoring. OPTIONAL PROTOCOL: Instruct the patient to void at 8 AM and discard the specimen. Then collect all urine including the final specimen voided at the end of the 24-hour collection period (ie, 8 AM the next morning). Avoid contact with metal during collection. Screw the lid on securely. See individual tests for collection times for industrial monitoring.  Special Instructions: Request form MUST state 24-hour specimen volume, if applicable.  Specimen Stability: Ambient: 14 days, Refrigerated: 14 days, Frozen: 14 days |                |   |                       |  |  |

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025

## Labcorp Oklahoma, Inc. Test Directory

| GENERAL INFORMATION             |  |                       |         |  |  |  |
|---------------------------------|--|-----------------------|---------|--|--|--|
| Expected TAT                    | 3 - 6 days   |                       |         |  |  |  |
| Performing Labcorp Test<br>Code | 007046   |                       |         |  |  |  |
| Notes                           | If total arsenic concentration is between 35-2000 ug/L, then Arsenic, Fractionated shown below, will be added to determine the proportion of organic, inorganic, and methylated forms. Additional charges apply. |                       |         |  |  |  |
|                                 | Arsenic, Fractionated [Order Code: 3709000]  |                       |         |  |  |  |
|                                 | RESULT CODE  | RESULT NAME           | LOINC   |  |  |  |
|                                 | 3800475  | Arsenic Inorganic     | 12481-8 |  |  |  |
|                                 | 3800276  | Arsenic, MMA          | 94458-7 |  |  |  |
|                                 | 3800476  | Arsenic, DMA          | 94457-9 |  |  |  |
|                                 | 3800676  | Arsenic, Total Toxic  | n/a     |  |  |  |
|                                 | 3800876  | Arsenic, Please Note: | n/a     |  |  |  |
|                                 | 3800875  | Arsenic Organic       | n/a     |  |  |  |
|                                 | 3800675  | Arsenic Methylated    | n/a     |  |  |  |
|                                 |  |                       |         |  |  |  |
| CPT Code(s)                     | 82570, 83825, 82175, 83655   |                       |         |  |  |  |
| Lab Section                     | Reference Lab  |                       |         |  |  |  |