

Trazodone Level

Order Name: **TRAZODONE**

Test Number: 3638625

Revision Date: 06/28/2023

TEST NAME	METHODOLOGY	LOINC CODE
Trazodone Level	Quantitative Liquid Chromatography/Tandem Mass Spectrometry	4064-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.6 mL)	Serum	Clot Activator (Red Top, No-Gel)	Room Temperature
Alternate 1	2 mL (0.6 mL)	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Room Temperature
Instructions	<p><b>Specimen Type:</b> Red-top tube OR green-top (heparin) tube. <b>DO NOT USE A GEL-BARRIER TUBE.</b> The use of gel-barrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant.</p> <p><b>Specimen Storage:</b> Room Temperature</p> <p><b>Specimen Collection:</b> Transfer separated serum or plasma to a plastic transport tube.</p> <p><b>Specimen Stability:</b> Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	2 - 4 days three working day.
Clinical Use	Trazodone is indicated for treatment of depression. Also known as: Desyrel.
Performing Labcorp Test Code	071688
Notes	Labcorp Test Code: 071688
CPT Code(s)	80338 (G0480)
Lab Section	Reference Lab