Labcorp Oklahoma, Inc. Test Directory

Tumor Necrosis Factor Alpha

Order Name: TNF

Test Number: 3657715
Revision Date: 12/12/2022

TEST NAME			METHODOLOGY	LOINC CODE
Tumor Necrosis Factor Alpha			Enzyme-Linked Immunosorbent Assay (ELISA)	3074-2
SPECIMEN REQU	UREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.4)	Serum	Clot Activator SST	Frozen
Alternate 1	1 mL (0.4)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	Notes: 0.4 mL (Note: This volume Does NOT allow for repeat testing.) Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: FREEZE Specimen Collection: Allow blood sample collected in serum separator tube to clot at room temperature for 30 minutes prior to centrifugation. Allowing samples to sit on the clot for more than 30 minutes may result in higher TNF-alpha levels. Transfer the serum into a LabCorp PP transpak frozen purple tube with screw cap LabCorp ID 23751). The specimen should be FROZEN immediately and maintained FROZEN until tested. To avoid delays in turnaround time when requesting multiple tests on frozen samples, please submit separate frozen specimens for each test requested.			

GENERAL INFORMATION		
Expected TAT	3 - 7 days	
Performing Labcorp Test Code	140673	
Notes	Labcorp Test Code: 140673	
CPT Code(s)	83520	
Lab Section	Reference Lab	

Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: 9 days

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025