Selenium Level Urine

Order Name: SELENIUM U Test Number: 3658675 Revision Date: 10/16/2024

TEST NAME			METHODOLOGY	LOINC CODE	
Creatinine Screen			Immunoassay (IA)	2161-8	
Selenium Level Uri	ne		Inductively-Coupled Plasma/Mass Spectrometry	5726-5	
Selenium/Creatinine		Calculation	13467-6		
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	

Preferred	3 mL (1.3 mL)	Urine, Random	Acid Washed, Trace Element Free Contatiner	Refrigerated		
Instructions	Specimen: 3mL(1.3mL) Random Urine in acid washed trace emelent free container					
	Collection: Unpreserved urine should be refrigerated immediately and analyzed within 1 week of collection. Acceptable preservatives include: Trace					
	Metal Free Hydrochloric Acid or Nitric Acid (0.1 mL of 12M acid/10 mL urine).					
	Stability Requirements: Room temperature n/a, Refrigerated 30 Days, Frozen 12 Months					

Cause for Rejection: Received Room Temperature

Testing referred to National Medical Services NMSPA#4180U

GENERAL INFORMATION								
Expected TAT	8 days (If Positive: 13 days) from set up							
Performing Labcorp Test Code	833870							
Notes								
	ORDER CODE	REFLEX NAME	RESULT CODE	RESULT NAME				
	4000073	.Rflex Specific Gravity	4000073	Specific Gravity Screen				
			1001500	Specific Gravity Urine				
CPT Code(s)	84255							
Lab Section	Reference Lab							