Testosterone Total, LCMS

Order Name: **TEST T LCMS** Test Number: 3602651 Revision Date: 06/28/2023

TEST NAME		METHODOLOGY	LOINC CODE	
Testosterone, Total	ne, Total		Liquid Chromatography/Tandem Mass S	pectrometry 2986-8
SPECIMEN REQUIRE	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator SST	Refrigerated
Instructions	 Recommended for Women, children under 13 years of age, hypogonadal males, or patients undergoing anti-androgenic therapies. Specimen: Serum - Red-top tube or serum separator (SST) Collection: Serum should be transferred from cells within one hour of collection and transferred to a plastic transport tube. Preferred Specimen Volume: 1.0 mL (0.2 mL min - does not allow for repeat testing) Storage/ Transport Temperature: Room temperature, refrigerated or frozen Stability: Room Temp: 8 days; Refrigerated: 7 days; Frozen: 5.58 years Cause for Rejection: Improperly labeled specimens, specimens collected tube other than red-top or SST. 			

GENERAL INFORMATION		
Testing Schedule	Mon-Fri	
Expected TAT	2-5 Days	
CPT Code(s)	84403	
Lab Section	Chemistry	