## Adenosine Deaminase, Pericardial Fluid

Order Name: Adeno D Pericard

Test Number: 6904899 Revision Date: 12/12/2022

TEST NAME		ا	METHODOLOGY	LOINC CODE	
Adenosine Deaminase, Pericardial Fluid			Quantitative Spectrophotometry	49760-2	
SPECIMEN REQU	IREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	0.5 mL (0.2 mL)	Pericardial Fluid	Sterile Screwtop Container	Frozen	
Instructions	Collect Pericardial Fluid, centrifuge specimen at room temperature. Transfer 0.5mL(0.1) fluid to a Standard Transport Tube and freeze. Indicate source on requisition.				
	Specimen must remain frozen until received in lab.				
	Unacceptable Conditions: Whole blood or Bronchoalveolar Lavage (BAL) specimens. Turbid specimens.  Special Instructions: Testing referred to ARUP TC 3002984				
	Stability: Ambient 24 hours; F	Refrigerated 7 days; Frozen 30 da	iys.		

GENERAL INFORMATION		
Expected TAT	8 - 13 days	
Performing Labcorp Test Code	830932	
Notes	Labcorp Test Code: 830932	
CPT Code(s)	84311	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2024