Labcorp Oklahoma, Inc. Test Directory

Zika Virus PCR, Blood

Order Name: Zika Virus PCR

Test Number: 6906615
Revision Date: 02/22/2018

TEST NAME			METHODOLOGY	LOINC CODE
Zika Virus PCR, Blood			Qualitative PCR	79190-5
SPECIMEN REQUIR	REMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1 mL)	Serum	Clot Activator SST	Frozen
Instructions	Separate from cells. Transfer 2 mL (Min: 1 mL) serum to a sterile container and Freeze. Unacceptable Conditions: Urine Remarks: Specimen source required. Stability: Ambient: Unacceptable; Refrigerated: 5 days; Frozen: 6 weeks			
	Please provide the following information on your requistion if possible Is the patient pregnant? - Relevent Zika Exposure? - Symptoms consistent with Zika?			
	Testing should only be performed on individuals meeting Centers for Disease Control and Prevention (CDC) clinical criteria for Zika virus (eg, clinical signs and symptoms associated with Zika virus infection) and/or CDC epidemiological criteria for Zika virus (eg, history of residence in or travel to a geographic region with active Zika transmission at the time of travel, or other epidemiological criteria for which Zika virus testing may be indicated).			

GENERAL INFORMATION		
Testing Schedule	Mon, Wed, Fri	
Expected TAT	1-4 days	
Clinical Use	Zika Virus by PCR is a real-time RT-PCR test intended for the qualitative detection of Zika virus RNA from individuals meeting CDC Zika virus clinical criteria (e.g., clinical signs and symptoms associated with Zika virus infection) and/or CDC Zika virus epidemiological criteria (e.g., history of residence in or travel to a geographic region with active Zika transmission at the time of travel, or other epidemiologic criteria for which Zika virus testing may be indicated).	
Notes	Prompt code: (6906619) Zika Virus B Source	
	Reference Lab: ARUP	
	Test Code: 2014065 Click Here to view information on the ARUP website.	
CPT Code(s)	87662	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2024