

## High-Sensitivity Troponin T

Order Name: **HS Troponin T**  
Test Number: 6907309  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
High-Sensitivity Troponin T	Quantitative Electrochemiluminescent Immunoassay	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.8 mL)	Plasma	Lithium Heparin PST (Light Green Top)	Frozen

**Instructions**

**Notes:** 0.8 mL (Note: This volume Does NOT allow for repeat testing).

**Specimen Type:** Green top (lithium-heparin) tube. Lithium-heparin plasma tubes containing separating gel can be used.

**Specimen Storage:** **FREEZE**

**Specimen Collection:** Centrifuge and transfer separated heparin plasma to a plastic transport tube. To avoid delays in turnaround time when requesting multiple tests on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.

**Special Instructions:** This test may exhibit interference when sample is collected from a person who is consuming a supplement with a high dose of biotin (also termed as vitamin B7 or B8, vitamin H, or coenzyme R). It is recommended to ask all patients who may be indicated for

**Specimen Stability:** Ambient: Not Available, Refrigerated : 24 hours (stability determined by manufacturer or, Frozen: 12 months (stability determined by manufacturer or

### GENERAL INFORMATION

Expected TAT	1 - 3 days
Performing Labcorp Test Code	140150
Notes	Labcorp Test Code: 140150
CPT Code(s)	84484
Lab Section	Reference Lab