

Folate, RBC

Order Name: **Folate RBC**
 Test Number: 5194359
 Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Folate, RBC		

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	See Instructions	See Instructions	See Instructions
Instructions	<p>Specimen Type: TWO full lavender-top (EDTA) tubes</p> <p>Specimen Storage: FREEZE Whole Blood in transport tube and Keep Room Temperature Whole Blood Specimen.</p> <p>Specimen Collection: Transfer 4 mL from One whole blood tube into a plastic transport tube and FREEZE. Keep the Second Tube Room temperature and whole blood - Do Not Spin.</p> <p>To avoid delays in turnaround time when requesting multiple tests on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST.</p> <p>Special Instructions: This test may exhibit interference when sample is collected from a person who is consuming a supplement with a high dose of biotin (also termed as vitamin B7 or B8, vitamin H, or coenzyme R). It is recommended to ask all patients who may be indicated for</p> <p>Specimen Stability: Ambient: Folate: Unstable Hematocrit: 1 day, Refrigerated : Folate: 1 day Hematocrit: 3 days, Frozen: Folate: 14 days Hematocrit: Unstable</p>			

GENERAL INFORMATION

Expected TAT	1 - 2 days
Performing Labcorp Test Code	266015
Notes	Labcorp Test Code: 266015
CPT Code(s)	85014, 82747
Lab Section	Reference Lab