

## Chlamydia/Gonococcus/Genital Mycoplasma Profile, NAA, Urine

Order Name: **Ct Ng Myco UR**  
 Test Number: 5195320  
 Revision Date: 09/12/2023

| TEST NAME   | METHODOLOGY                      | LOINC CODE |
|---|----------------------------------|------------|
| Chlamydia/Gonococcus/Genital Mycoplasma Profile, NAA, Urine | Nucleic acid amplification (NAA) |            |

### SPECIMEN REQUIREMENTS

| Specimen  | Specimen Volume (min) | Specimen Type     | Specimen Container          | Transport Environment |
|-----------|-----------------------|-------------------|-----------------------------|-----------------------|
| Preferred | 30 mL (5 mL)          | Urine (See Inst.) | APTIMA Urine Transport Tube | Room Temperature      |

**Instructions**

**Specimen:** First-void urine in Aptima® urine transport  
**Volume:** 2 mL of a 20 to 30 mL urine collection  
**Container:** Aptima® urine transport  
**Collection:** Urine specimen: The patient should not have urinated for at least one hour prior to specimen collection. Direct patient to provide a first-catch urine (approximately 20 to 30 mL of the initial urine stream) into a urine collection cup free of any preservatives. Collection of larger volumes of urine may result in specimen dilution that may reduce test sensitivity; lesser volumes may not adequately rinse organisms into the specimen. Female patients should not cleanse the labial area prior to providing the specimen. Add urine to the Aptima® Combo 2 urine collection device. The final volume must be between the two black lines on the device (about 2 mL).  
**Storage Instructions:** Maintain specimen at room temperature or refrigerate (2°C to 30°C).  
**Stability** Room temperature 30 days (stability provided by manufacturer or literature reference)

### GENERAL INFORMATION

|                              |          |
|------------------------------|----------|
| Expected TAT                 | 3-4 days |
| Performing Labcorp Test Code | 180051   |

| Notes   | RESULT CODE      | RESULT NAME    | LOINC   |
|---------|------------------|----------------|---------|
|         | 6907667          | M genitalia Ur | 88226-6 |
| 5195379 | M. hominis Urine | 68546-1        |         |
| 5195380 | Ureaplasma Ur    | 68547-9        |         |
| 6910151 | Chlamydia Ur     | 42931-6        |         |
| 6910155 | Gonorrhea Ur     | 60256-5        |         |

|             |                              |
|-------------|------------------------------|
| CPT Code(s) | 87491, 87591, 87798x2, 87563 |
| Lab Section | Reference Lab                |