Vitamin B6

Order Name: VIT B6 Test Number: 3603660 Revision Date: 12/12/2022

TEST NAME			METHODOLOGY	LOINC CODE	
Vitamin B6			Liquid Chromatography/Tandem Mass S	Spectrometry 62236-5	
SPECIMEN REQUIREM	ENTS				
		Specimen Tune	Specimen Container	Transport Environment	
Specimen	Specimen Volume (min)	Specimen Type	· · · · · · · · · · · · · · · · · · ·		
Preferred	1 mL (0.25)	Plasma	EDTA (Lavender Top)	Frozen (Light Protect)	
Instructions	tube completely, top and bottom, aluminum foil. Secure with tape.) Container Detail: Protect from L Specimen Storage: Refrigerate Specimen Collection: Collect b six times to ensure adequate mix Specimens should be stored refr delays in turnaround time when r TEST REQUESTED.	 cimen Type: Lavender-top (EDTA) tube; amber plastic transport tube with amber-top. (If amber tubes are unavailable, cover standard transport completely, top and bottom, with aluminum foil. Identify specimen with patient's name directly on the container AND on the outside of the inum foil. Secure with tape.) For amber plastic transport tube and amber-top, order LabCorp item No. 23594. ainer Detail: Protect from Light cimen Storage: Refrigerated or FROZEN and PROTECT FROM LIGHT. cimen Collection: Collect blood by venipuncture into a lavender-top tube containing EDTA and mixed immediately by gentle inversion at least mes to ensure adequate mixing. The plasma must be separated and protected from light in an amber transport tube with amber stopper. imens should be stored refrigerated or frozen immediately and maintained at temperature during shipping and at the testing facility. To avoid rs in turnaround time when requesting multiple tests on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH 			
GENERAL INFORMATION					
Expected TAT 4 - 6 days week day.					
Clinical Use	Vitamin B6 is a cofactor in many metabolic pathways including heme synthesis. Vitamin B6 deficiency may be observed in patients with metabolic disorders, secondary to therapeutic drug use, or alcoholism. Deficiency affects the function of the immune system.				
Performing Labcorp Tes Code	t 004655				
Notes	Labcorp Test Code: 004	4655			
CPT Code(s)	84207	84207			
Lab Section	Reference Lab				