

Deoxycorticosterone

Order Name: **DEOXYSTER**

Test Number: 3603055

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Deoxycorticosterone	Mass Spectrometry (MS)	1656-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3mL (0.5)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	<p><b>Notes:</b> 0.5 mL (Note: This volume Does NOT allow for repeat testing.)</p> <p><b>Specimen Type:</b> Red-top tube, gel-barrier tube, lavender-top (EDTA) tube OR green-top (heparin) tube</p> <p><b>Specimen Storage:</b> <b>FREEZE</b></p> <p><b>Specimen Collection:</b> <b>Separate serum from cells within one hour of collection. Transfer specimen to a plastic transport tube and Freeze.</b> To avoid delays in turnaround time when requesting multiple test on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.</p> <p><b>Special Instructions:</b> Testing is performed at Esoterix Endocrinology Laboratory: ES#500124</p> <p><b>Specimen Stability:</b> Ambient: 3 days, Refrigerated : 3 days, Frozen: 2 years</p>			

GENERAL INFORMATION	
Expected TAT	6 - 15 days
Performing Labcorp Test Code	500138
Notes	Labcorp Test Code: 500138
CPT Code(s)	82633
Lab Section	Reference Lab