## Labcorp Oklahoma, Inc. Test Directory

## Pneumococcal Antibody Panel (23 Serotype)

Order Name: PNEUMO 23
Test Number: 5575605
Revision Date: 01/01/2025

TEST NAME			METHODOLOGY	LOINC CODE	
Pneumococcal Antibody Panel (23 Serotype)			Fluoroimmunoassay		
SPECIMEN REQU	UREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	1.0 mL (0.5 mL)	Serum	Clot Activator SST	Room Temperature	
Instructions	Notes: 0.5 mL (Note: This volume Does NOT allow for repeat testing).				
	Specimen Type: Gel-barrier tube				
	Specimen Storage: Maintain specimen at room temperature.				
	Specimen Collection: Not Available				
	Special Instructions: Testing referred to Eurofins Viracor LLC EURKS#401835				

Specimen Stability: Ambient: 1 week, Refrigerated: 1 week, Frozen: 5 months

GENERAL INFORMATION				
Expected TAT	7-10 Days			
Clinical Use	Includes Serotypes 1, 2, 3, 4, 5, 6B, 7F, 8, 9N, 9V, 10A, 11A, 12F, 14, 15B, 17F, 18C, 19A, 19F, 20, 22F, 23F, and 33F. (Conjugated Serotypes: 1, 3, 4, 5, 6B, 7F, 9V, 14, 18C, 19A, 19F, and 23F.).			
	A pre-and post-vaccination sample comparison is required in order to assess the humoral immune response to vaccination with Streptococcus pneumoniae vaccine. Post vaccination samples should be obtained at 4 weeks post-immunization. A single sample provides only general immune status of an individual to various serotypes of pneumococci. Immune response to pneumococci may vary in different individuals based on the age, vaccination status, immunologic competence and the serotype of the organism. Long-term protection is generally considered to be a serotype antibody level of 1.3 ug/ mL or greater.			
	Following vaccination, an antibody ratio of less than twofold is considered a non-responder; a ratio of two-to fourfold is a weak responder; a ratio of fourfold or greater is a good responder. The higher the pre-vaccination antibody level for a specific pneumococcal serotype, the less likely the response will increase significantly after vaccination.			
	An increased antibody level to 50-70 percent or more of the serotypes is thought to represent a normal humoral response. In the case of pure polysaccharide vaccine, indication of immune system competence is further delineated as an adequate response to at least 70 percent of the serotypes in the vaccine challenge for those 6-65 years of age, or to at least 50 percent of the serotypes in the vaccine challenge for those 2-5 years of age.			
Performing Labcorp Test Code	812166			
CPT Code(s)	86581			
Lab Section	Reference Lab			

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025