Labcorp Oklahoma, Inc. Test Directory

Toluene, Blood

Order Name: **TOLUENE**Test Number: 3610275
Revision Date: 12/12/2022

LOINC

TEST NAME			METHODOLOGY	CODE	
Toluene, Blood		Headspace Gas Chromatography/Flame Ionizati FID)	ion Detection (GC- 9637-0		
SPECIMEN REQUIR	REMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	20 mL (7 mL)	Whole Blood	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated	
Instructions	Specimen Type: Green-top (heparin) tube				
	Specimen Storage: Refrigerate; specimen may also be frozen.				
	Specimen Collection: Collect adequate number of tubes to obtain the required specimen volume. Tubes should be completely filled to prevent loss				
	of volatile compound into headspace. Collect prior to last shift of workweek.				
	Special Instructions: Testing referred to MedTox Laboratories Inc TC 70123				

GENERAL INFORMATION		
Expected TAT	10 - 14 days	
Performing Labcorp Test Code	790580	
Notes	Labcorp Test Code: 790580	
CPT Code(s)	84600	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025