Pyruvic Acid, CSF

Order Name: PYRU CSF
Test Number: 6905131
Revision Date: 12/12/2022

TEST NAME		METI	HODOLOGY	LOINC CODE
Pyruvic Acid, CSF		Quantitative Enzymatic		14122-6
SPECIMEN REQUIRE	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	CSF (Cerebrospinal Fluid)	Sterile Screwtop Container	Frozen
Instructions	Specimen Type: Sterile screw capped vial Specimen Storage: Freeze as soon as possibe after collection Specimen Collection: Collect entire sample into a single sterile tube. Immediately after collection, add a 2:1 ratio (acid:sample) of COLD 8% perchloric acid. Vortex for 30 seconds, and place on ice for 5 minutes. Store frozen at -20 degrees Celsius and ship frozen. Special Instructions: Testing referred to Medical Neurogenetic LLC MNEGA#MET11 Specimen Stability: Ambient: unstable (stability provided by manufacturer or li, Refrigerated : 24 hours (stability provided by manufacturer or li, Frozen: -20 C = 72 hours; -80 C = Indefinate (stability pr			

GENERAL INFORMATION		
Expected TAT	10 - 14 days	
Performing Labcorp Test Code	620045	
Notes	Labcorp Test Code: 620045	
CPT Code(s)	84210	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025