

Sertraline, Serum or Plasma

Order Name: Sertraline

Test Number: 7001007

Revision Date: 11/14/2022

TEST NAME	METHODOLOGY	LOINC CODE
Sertraline, Serum or Plasma	Quantitative Liquid Chromatography/Tandem Mass Spectrometry	6906-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	<p>Specimen Type: Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.</p> <p>Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p>Specimen Collection: Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.</p> <p>Special Instructions: Testing referred to MedTox Laboratories Inc MX#122</p>			

GENERAL INFORMATION	
Expected TAT	5-10 days
Performing Labcorp Test Code	808494
Notes	Labcorp Test Code: 808494
CPT Code(s)	80299
Lab Section	Reference Lab