

High-Sensitivity Troponin T

Order Name: **HS Troponin T**

Test Number: 6907309

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
High-Sensitivity Troponin T	Quantitative Electrochemiluminescent Immunoassay	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.8 mL)	Plasma	Lithium Heparin PST (Light Green Top)	Frozen
Instructions	<p><b>Notes:</b> 0.8 mL (Note: This volume Does NOT allow for repeat testing).</p> <p><b>Specimen Type:</b> Green top (lithium-heparin) tube. Lithium-heparin plasma tubes containing separating gel can be used.</p> <p><b>Specimen Storage:</b> <b>FREEZE</b></p> <p><b>Specimen Collection:</b> Centrifuge and transfer separated heparin plasma to a plastic transport tube. To avoid delays in turnaround time when requesting multiple tests on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.</p> <p><b>Special Instructions:</b> This test may exhibit interference when sample is collected from a person who is consuming a supplement with a high dose of biotin (also termed as vitamin B7 or B8, vitamin H, or coenzyme R). It is recommended to ask all patients who may be indicated for</p> <p><b>Specimen Stability:</b> Ambient: Not Available, Refrigerated : 24 hours (stability determined by manufacturer or, Frozen: 12 months (stability determined by manufacturer or</p>			

GENERAL INFORMATION	
Expected TAT	1 - 3 days
Performing Labcorp Test Code	140150
Notes	Labcorp Test Code: 140150
CPT Code(s)	84484
Lab Section	Reference Lab