Arsenic 24 Hour Urine or Random Urine

Order Name: ARSENIC U
Test Number: 3809050
Revision Date: 03/27/2025

TEST NAME			METHODOLOGY	LOINC CODE		
Arsenic 24 Hour Urine or Random Urine			Inductively Coupled Plasma/Mass Spectrometry (ICP/MS) or Atomic Spectroscopy (AS)			
SPECIMEN REQUI	REMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment		
Preferred	5 mL (2.5 mL)	Urine, 24-hour	24 hour Urine Container	Room Temperature		
Alternate 1	5 mL (2.5 mL)	Urine, Random	Sterile Urine container	Room Temperature		
Instructions	Specimen: Urine (random or 24-hour) Notes: 2.5 mL (Note: This volume Does NOT allow for repeat testing.) Specimen Type: Plastic urine container, no preservative Specimen Storage: Maintain specimen at room temperature. Specimen Collection: Sampling time is the end of the work week for industrial exposure monitoring. Metals with timing "end of work week" (meaning four or five consecutive working days with exposure) are eliminated with half-lives longer than five hours. Such metals accumulate in the body during the work week; therefore, their timing is critical in relation to previous exposures. OPTIONAL PROTOCOL: Instruct the patient to void at 8 AM and discard the specimen. Then collect all urine including the final specimen voided at the end of the 24-hour collection period (ie, 8 AM the next morning). Avoid contact with metal during collection. Screw the lid on securely. Special Instructions: Request form MUST state volume collected as 24-hour, if applicable. Do not use preservatives. Preservatives used for routine analysis may contain mercuric oxide (ie, Stabilur), which interferes with all metal testing. If both urinalysis and metal testing					

Specimen Stability: Ambient: 14 days, Refrigerated: 14 days, Frozen: 14 days (Freeze/thaw cycles Stable x3)

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025

Labcorp Oklahoma, Inc. Test Directory

GENERAL INFORMATION						
Expected TAT	3 - 6 days					
Performing Labcorp Test Code	007045					
Notes	If total arsenic concentration is between 35-2000 ug/L, then Arsenic, Fractionated shown below, will be added to determine the proportion of organic, inorganic, and methylated forms. Additional charges apply.					
	Arsenic, Fractionated [Order Code: 3709000]					
	RESULT CODE	RESULT NAME	LOINC			
	3800475	Arsenic Inorganic	12481-8			
	3800276	Arsenic, MMA	94458-7			
	3800476	Arsenic, DMA	94457-9			
	3800676	Arsenic, Total Toxic	n/a			
	3800876	Ars(Toxic)/Crt Note:	n/a			
	3800875	Arsenic Organic	n/a			
	3800675	Arsenic Methylated	n/a			
CPT Code(s)	82570, 82175					
Lab Section	Reference Lab					