

Lupus Anticoagulant-Hexagonal Phospholipid Neutralization

Order Name: **HEXA PHOS**

Test Number: 1507375

Revision Date: 06/21/2018

TEST NAME	METHODOLOGY	LOINC CODE
Lupus Anticoagulant-Hexagonal Phospholipid Neutralization	Clot Detection	33930-9

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5.4 mL	Whole Blood	Sodium Citrate 3.2% (Blue Top)	Room Temperature
Alternate 1	3.0 mL	Double Spun Plasma	Sterile, Capped Plastic Tube	Frozen
Instructions	Please indicate anticoagulant therapy. Collect Two 2.7mL Sodium Citrate 3.2% (Blue Top) tubes. Each 2.7mL Sodium Citrate 3.2% (Blue Top) tube must be filled to the proper level, no hemolysis. Improperly filled tubes can give erroneous results. Whole blood must be transported to lab immediately. If testing cannot be started within 4 hours of collection the specimen must be double spun then 1.5mL plasma aliquot from each tube into individual plastic aliquot tubes and freeze. Do not pool aliquots together!			

GENERAL INFORMATION	
Testing Schedule	Tues, Thur
Expected TAT	7-10 Days
CPT Code(s)	85598
Lab Section	Coagulation