

Smooth Muscle Antibody Titer

Order Name: **SM MUSC DE**

Test Number: 5572400

Revision Date: 10/23/2017

TEST NAME	METHODOLOGY	LOINC CODE
Smooth Muscle Antibody Titer	Indirect Fluorescent Antibody	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.0 mL	Serum	Clot Activator SST	Refrigerated
Instructions	Allow specimen to clot completely at room temperature. Separate serum or plasma from cells ASAP or within 2 hours of collection. Stability After separation from cells: Ambient 4hours, Refrigerated 7 days, Frozen 1 month (avoid repeated freeze/thaw cycles).			

GENERAL INFORMATION

Testing Schedule	Mon - Fri
Expected TAT	3 Days
Clinical Use	Smooth muscle antibody is in high titers (\geq 1:160) in approximately 97% of patients with autoimmune chronic active hepatitis.
CPT Code(s)	86256
Lab Section	Immunology - Serology