

## Cyclosporine Level Peak

Order Name: **CYCLO PEAK**

Test Number: 2015325

Revision Date: 05/03/2021

TEST NAME	METHODOLOGY	LOINC CODE
Cyclosporine Level Peak	Chemiluminescent Micropartical Immunoassay	53834-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>2 mL (1.0)</b>	<b>Whole Blood</b>	<b>EDTA (Lavender Top)</b>	<b>Refrigerated</b>
Instructions	Peak: draw 2 hours after medication. Stability: Ambient 24 hours. Refrigerated 7 days. Frozen > 7 days.			

GENERAL INFORMATION	
Testing Schedule	Mon - Fri
Expected TAT	1-2 days
Clinical Use	This test is used to monitor blood concentration of Cyclosporine A.
CPT Code(s)	80158
Lab Section	Chemistry