

## Urea Nitrogen Urine Random

Order Name: **UUN R U**  
Test Number: 3000270  
Revision Date: 06/04/2021

TEST NAME	METHODOLOGY	LOINC CODE
Urea Nitrogen Urine Random	Urease/GLDH	3095-7

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Random urine collection. No preservative. Keep refrigerated. Specimen stability: Ambient 2 days, Refrigerated 7 days, Frozen 1 month			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful as an aid in evaluating renal function.
CPT Code(s)	84540
Lab Section	Chemistry