Labcorp Oklahoma, Inc. Test Directory

Vancomycin Peak

Order Name: VANCO PEAK

Test Number: 4005900 Revision Date: 06/04/2021

TEST NAME			METHODOLOGY	LOINC CODE
Vancomycin Peak			Enzyme Multiplied Immunoassay Technique	4090-7
SPECIMEN REQUI	REMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated or Frozen
Alternate 1	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated or Frozen
Instructions	Separate plasma/serum from o	Peak: draw specimen 60 minutes after a 1 hour infusion. Separate plasma/serum from cells/gel immediately. Stability: Refrigerated 7 days, Frozen 14 days.		

GENERAL INFORMATION		
Testing Schedule	Daily	
Expected TAT	1-2 days	
Clinical Use	Useful for optimizing drug dosage and assessing toxicity.	
CPT Code(s)	80202	
Lab Section	Chemistry	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025