Labcorp Oklahoma, Inc. Test Directory

RBC Osmotic Fragility

Order Name: FRAGILITY
Test Number: 0104740
Revision Date: 12/12/2022

TEST NAME		METH	ODOLOGY	LOINC CODE
RBC Osmotic Fragility		Spect	pectrophotometry	
SPECIMEN REQUIREM	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5mL + 2 Smears (1mL + 2 Smears)	Whole Blood	Sodium Heparin (Green Top / No- Gel)	Refrigerated
Alternate 1	5mL + 2 Smears (1mL + 2 Smears)	Whole Blood	EDTA (Lavender Top)	Refrigerated
Alternate 2	5mL (1mL + Smears)	Whole Blood	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated
Instructions	Notes: 1 mL (Note: This volume Does NOT allow for repeat testing). Specimen Type: Green-top (heparin) tube OR lavender-top (EDTA) tube AND (2) unfixed smears Specimen Storage: Refrigerate Stability: Room temperature: unacceptable Refrigerate: 72 hours Frozen: unacceptable Specimen Collection: Draw specimen Mon-Wed ONLY. Refrigerate specimens within 30 minutes of collection. Specimens must be received and testing performed by ARUP within 72 hours of collection. Special Instructions: Testing referred to ARUP Y8#2002257 NOTE: Draw specimen Mon-Wed ONLY. Please contact regional referral department to arrange transportation of time sensitive specimen Specimen Stability: Ambient: Not Available, Refrigerated: Not Available, Frozen: Not Available			

GENERAL INFORMATION		
Expected TAT	5 - 8 days	
Clinical Use	To confirm suspected red cell spherocytosis.	
Performing Labcorp Test Code	284205	
Notes	Labcorp Test Code: 284205	
CPT Code(s)	85555; 80500	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025