

Prostatic Acid Phosphatase (PAP)

Order Name: **ACID PHS P**

Test Number: 2000275

Revision Date: 04/28/2023

TEST NAME	METHODOLOGY	LOINC CODE
Prostatic Acid Phosphatase (PAP)	Immunochemiluminometric (ICMA)	20420-6

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.5 mL (0.2)	Serum	Clot Activator (SST or Red No-Gel)	Frozen
Instructions	<p>Notes: 0.2 mL (Note: This volume Does NOT allow for repeat testing.)</p> <p>Specimen Type: Red-top tube or gel-barrier tube; do NOT use acid phosphatase transport tube.</p> <p>Specimen Collection: Transfer the serum into a Labcorp PP transpak frozen purple tube with screw cap (Labcorp No.49482) or Sterile Screwtop container. Freeze immediately and maintain frozen until tested. To avoid delays in turnaround time when requesting multiple tests on frozen samples, please submit separate frozen specimens for each test requested.</p> <p>Causes for Rejection: Plasma specimen; specimen not frozen.</p> <p>Specimen Stability: Ambient: 3 hrs; Refrigerated: 24 hrs; Frozen: 6 months</p>			

GENERAL INFORMATION	
Expected TAT	3-4 days
Clinical Use	This test is not the same as Acid Phosphatase. Values obtained with different assays should not be used interchangeably in serial testing. It is recommended that only one assay method be used consistently to monitor each patient's course of therapy. This procedure does not provide serial monitoring; it is intended for one-time use only.
Performing Labcorp Test Code	004747
Notes	Labcorp Test Code: 004747
CPT Code(s)	84066
Lab Section	Reference Lab