

Latex Allergen

Order Name: **LATEX IGE**

Test Number: 5610780

Revision Date: 10/23/2017

TEST NAME	METHODOLOGY	LOINC CODE
Latex Allergen	ImmunoCAP	15806-3

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.4 mL (0.2 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	<p>Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container.</p> <p>Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)</p>			

GENERAL INFORMATION	
Expected TAT	3-5 days
Clinical Use	This test can be used to detect a general allergy to all latex products.
Notes	Tests for Latex Hevea brasiliensis (K82) the Non-ammoniated, Buffered Latex
CPT Code(s)	86003
Lab Section	Reference Lab