Labcorp Oklahoma, Inc. Test Directory

Beta-2, Transferrin

Preferred

Instructions

Order Name: **BETA-2 TRA**Test Number: 3656675
Revision Date: 12/12/2022

Refrigerated

Sterile Screwtop Container

| TEST NAME | | | METHODOLOGY | LOINC CODE |
|-------------------|-----------------------|---------------|--------------------------------|-----------------------|
| Beta-2, Transferr | in | | Immunofixation Electrophoresis | |
| SPECIMEN REQ | UIREMENTS | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |

Aural or nasal fluid in sterile container; Do Not Freeze. Do Not send CSF specimens!

Notes: 1 mL (Note: This volume Does NOT allow for repeat testing).

Specimen Type: Sterile container **Specimen Storage:** Refrigerated

2mL (1mL)

Special Instructions: Testing referred to ARUP TC 50047

Specimen Stability: Ambient: 2 weeks; Refrigerated: 2 weeks; Frozen: 2 weeks

Fluid

| GENERAL INFORMATION | | |
|---------------------------------|--|--|
| Expected TAT | 4 - 8 days | |
| Clinical Use | This test is to detect CSF in body fluids such as sinus or ear drainage. | |
| Performing Labcorp Test Code | 829030 | |
| Notes | Labcorp Test Code: 829030 | |
| CPT Code(s) | 86335 | |
| Lab Section | Reference Lab | |
| | | |

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025