Labcorp Oklahoma, Inc. Test Directory

Beta-2, Transferrin

Preferred

Instructions

Order Name: **BETA-2 TRA**Test Number: 3656675
Revision Date: 12/12/2022

Refrigerated

Sterile Screwtop Container

		METHODOLOGY	LOINC CODE
		Immunofixation Electrophoresis	
ENTS			
Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
		ENTS	Immunofixation Electrophoresis

Aural or nasal fluid in sterile container; Do Not Freeze. Do Not send CSF specimens!

Notes: 1 mL (Note: This volume Does NOT allow for repeat testing).

Fluid

Specimen Type: Sterile container **Specimen Storage:** Refrigerated

2mL (1mL)

Special Instructions: Testing referred to ARUP TC 50047

Specimen Stability: Ambient: 2 weeks; Refrigerated: 2 weeks; Frozen: 2 weeks

GENERAL INFORMATION		
Expected TAT	4 - 8 days	
Clinical Use	This test is to detect CSF in body fluids such as sinus or ear drainage.	
Performing Labcorp Test Code	829030	
Notes	Labcorp Test Code: 829030	
CPT Code(s)	86335	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2024