

## Sirolimus (Rapamycin)

Order Name: **SIROLIMUS**

Test Number: 3658510

Revision Date: 06/04/2021

| TEST NAME             | METHODOLOGY             | LOINC CODE |
|-----------------------|-------------------------|------------|
| Sirolimus (Rapamycin) | Chemiluminescence Assay | 29247-4    |

| SPECIMEN REQUIREMENTS |   |               |                     |                       |
|-----------------------|---|---------------|---------------------|-----------------------|
| Specimen              | Specimen Volume (min)   | Specimen Type | Specimen Container  | Transport Environment |
| Preferred             | 3 mL (2)  | Whole Blood   | EDTA (Lavender Top) | Refrigerated          |
| Instructions          | Ship and store refrigerated. Shipping at ambient temperature (< 25 C) is acceptable.<br>Optimal time to draw specimen: 1/2 to 1 hour before next oral dose at steady state.<br>Stability: Ambient n/a, Refrigerated 7 days. Frozen > 7 days |               |                     |                       |

| GENERAL INFORMATION |                             |
|---------------------|-----------------------------|
| Testing Schedule    | Mon - Fri                   |
| Expected TAT        | 1-3 Days                    |
| Notes               | *{ Note: 2006 CPT Updated.} |
| CPT Code(s)         | 80195                       |
| Lab Section         | Chemistry                   |