

## Estriol, Serum

Order Name: **ESTRIOL**  
Test Number: 3801500  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Estriol, Serum	Immunochemiluminometric (ICMA)	2251-7

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.3)	Serum	Clot Activator SST	Refrigerated
<b>Instructions</b>	<p><b>Notes:</b> 0.3 mL (Note: This volume Does NOT allow for repeat testing.)  <b>Specimen Type:</b> Red-top tube or gel-barrier tube  <b>Specimen Storage:</b> Refrigerate  <b>Specimen Collection:</b> If red-top tube is used, transfer separated serum to a plastic transport tube. All specimens should be drawn at the same time of day to allow comparison of values.  <b>Special Instructions:</b> State weeks of gestation on the request form. Usually done as serial measurement.  <b>Specimen Stability:</b> Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p>			

### GENERAL INFORMATION

Expected TAT	2 - 4 days
Performing Labcorp Test Code	004614
Notes	Labcorp Test Code: 004614
CPT Code(s)	82677
Lab Section	Reference Lab