

## Lorazepam (Ativan)

Order Name: **LORAZEPAM**

Test Number: 3611975

Revision Date: 12/20/2022

TEST NAME	METHODOLOGY	LOINC CODE
Lorazepam (Ativan)	Quantitative Liquid Chromatography/Tandem Mass Spectrometry	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>3 mL (0.6 mL)</b>	<b>Serum</b>	<b>Clot Activator SST</b>	<b>Room Temperature</b>
Instructions	<p><b>Specimen Type:</b> Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.</p> <p><b>Specimen Storage:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p><b>Specimen Collection:</b> Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.</p> <p><b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc MX#64</p>			

### GENERAL INFORMATION

Expected TAT	7-10 Days
Performing Labcorp Test Code	808465
Notes	Labcorp Test Code: 808465
CPT Code(s)	80299
Lab Section	Reference Lab