

## Helicobacter pylori Antigen Detection Stool

Order Name: **HPYLORI AG**

Test Number: 3502325

Revision Date: 09/21/2022

| TEST NAME                                   | METHODOLOGY      | LOINC CODE |
|---|------------------|------------|
| Helicobacter pylori Antigen Detection Stool | Immunoassay (IA) | 80373-4    |

| SPECIMEN REQUIREMENTS |   |                      |                                 |                       |
|-----------------------|---|----------------------|---------------------------------|-----------------------|
| Specimen              | Specimen Volume (min)   | Specimen Type        | Specimen Container              | Transport Environment |
| Preferred             | <b>5g (1)</b>   | <b>Stool, Random</b> | <b>Stool specimen container</b> | <b>Refrigerated</b>   |
| <b>Instructions</b>   | <b>Patient Preparation:</b> To avoid false negative results, no antibiotics, proton pump inhibitors, or a bismuth preparation should be administered to patient for 14 days prior to specimen collection.<br><br>Transfer 5g(1g) stool to an unpreserved stool transport vial and Freeze.<br><b>Unacceptable Conditions:</b> Gastric specimens, swabs or tissue. Specimens in media or preservatives.<br><b>Stability:</b> Ambient: 2 hours; Refrigerated: 72 hours; Frozen: 1 month. |                      |                                 |                       |

| GENERAL INFORMATION |              |
|---------------------|--------------|
| Testing Schedule    | Sun-Sat      |
| Expected TAT        | 2 Days       |
| CPT Code(s)         | 87338        |
| Lab Section         | Microbiology |