

# Cryofibrinogen

Order Name: **CRYOFIBRIN**  
 Test Number: 5221675  
 Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Cryofibrinogen	Precipitation	

## SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3mL (1mL)	Plasma	Sodium Citrate 3.2% (Blue Top)	See Instructions
Instructions	<p><b>Specimen Type:</b> Blue-top (sodium citrate) tube, lavender-top (EDTA) tube, or gray-top (sodium fluoride) tube  <b>Specimen Storage:</b> Maintain specimen at room temperature.  <b>Specimen Collection:</b> Specimen must be drawn in a prewarmed tube and kept at 37(degrees)C. Separate plasma from cells immediately and transfer the plasma into a clean transport tube. LABEL TRANSPORT TUBE "PREWARMED."  <b>Specimen Stability:</b> Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p>			

## GENERAL INFORMATION

Expected TAT	3 - 5 days later.
Performing Labcorp Test Code	080317
Notes	Labcorp Test Code: 080317
CPT Code(s)	82585
Lab Section	Reference Lab