

## Wet Prep for Yeast and Trichomonas

Order Name: **Wet Prep**  
Test Number: 2915475  
Revision Date: 11/15/2021

TEST NAME	METHODOLOGY	LOINC CODE
Vaginal Yeast Examination	Microscopy	32765-0
Trichomonas Microscopic Examination	Microscopy	32766-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	See Instructions	Wet prep saline	Room Temperature
Alternate 1	Swab	Swab	Copan eSwab - White (Regular size)	Room Temperature
Instructions	Collect discharge and place in sterile tube with saline or Copan eSwab - White. Deliver specimens to the laboratory for testing as soon as possible, ideally within 2 hours of collection. Trichomonas viability decreases in the specimen overtime. The <b>Trichomonas antigen (2915460)</b> testing is strongly recommended. Specimens arriving 24 hours after collection are unacceptable.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1 Day
Clinical Use	Determines yeast vaginitis and/or trichomoniasis
Notes	This used to be known as: Order Name <b>C WET PR</b> Test Number <b>6000650</b>
CPT Code(s)	87210
Lab Section	Microbiology