Alcohol Plasma/Serum

Order Name: ALCOHOL Test Number: 4000800 Revision Date: 04/26/2021

TEST NAME		м	IETHODOLOGY		
Alcohol Plasma/Serum		E	Enzymatic	5643-2	
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	1 mL (0.1 mL)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated	
Alternate 1	1 mL (0.1 mL)	Serum	Clot Activator SST	Refrigerated	
Instructions	Use betadine or phisohex prep sponge. Do not use alcohol. MUST be submitted in unopened original tube as exposure to air will affect result. Stability: Ambient 2 weeks, Refrigerated 6 months, Frozen 6 months				

GENERAL INFORMATION		
Testing Schedule	Daily	
Expected TAT	24hrs	
CPT Code(s)	G0480	
Lab Section	Chemistry	