## Oak Tree Allergen

Order Name: OAK TREE
Test Number: 5604925
Revision Date: 12/05/2016

TEST NAME			METHODOLOGY	LOINC CODE
Oak Tree Allergen			ImmunoCAP	6189-5
SPECIMEN REQUIRE	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.4 mL (0.2 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)			

GENERAL INFORMATION		
Expected TAT	3-5 days	
CPT Code(s)	86003	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025