## Bakers Yeast Allergen IgG

Order Name: YEAST IGG Test Number: 3666150 Revision Date: 12/14/2022

TEST NAME			METHODOLOGY	LOINC CODE	
Bakers Yeast Allergen IgG			ImmunoCAP	35538-8	
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	0.3 mL (0.1 mL)	Serum	Clot Activator SST	Room Temperature	
Instructions	Specimen Type: Red-top tube or gel-barrier tube				
	Specimen Storage: Room ter	nperature			
	Specimen Collection: Not Av	ailable			
	Special Instructions: Values obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only				
	one assay method be used consistently to monitor each patient's course of therapy.				
	Specimen Stability: Ambient:	14 days, Refrigerated : 14 day	s, Frozen: 3 months		
GENERAL INFORMATION					

Expected TAT	2 - 4 days	
Performing Labcorp Test Code	605066	
Notes	Labcorp Test Code: 605066	
CPT Code(s)	86001	
Lab Section	Reference Lab	