

## Burkitt's Lymphoma/NHL/ALL, IGH/MYC, t(8;14) by FISH

Order Name: **BURKI FSH**

Test Number: 9117900

Revision Date: 10/01/2022

TEST NAME	METHODOLOGY	LOINC CODE
Burkitt's Lymphoma/NHL/ALL, IGH/MYC, t(8;14) by FISH	Fluorescence in Situ Hybridization	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL (1 mL)	Bone Marrow	Sodium Heparin (Green Top, No-Gel)	Room Temperature
Alternate 1	5 mL (3 mL)	Whole Blood	Sodium Heparin (Green Top, No-Gel)	Room Temperature
Instructions	Send specimen ASAP, Keep at room temperature! (DO NOT FREEZE). Frozen samples will be rejected.			

### GENERAL INFORMATION

Testing Schedule	Mon-Fri
Expected TAT	3-5 Days
Clinical Use	Useful to diagnose Burkitt's type lymphoma and some cases of acute lymphoblastic leukemia.
CPT Code(s)	88377
Lab Section	Reference Lab