

Lambert-Eaton Myasthenic Syndrome (LEMS) Panel

Order Name: **LAMB EATON**

Test Number: 5503127

Revision Date: 06/01/2015

TEST NAME	METHODOLOGY	LOINC CODE
Striated Muscle Antibody Screen with Titer	Immunofluorescence Assay (IFA)	5372-8
Voltage-Gated Calcium Channel (VGCC) Antibody Assay	Quantitative Radioimmunoassay	31024-3
Acetylcholine Receptor Binding Antibody	See Test Notes	11034-6
Acetylcholine Receptor Modulating Antibody	Flow cytometry	30192-9

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	4 mL (2 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	Separate serum from cells ASAP or within 2 hours of collection.			

GENERAL INFORMATION

Testing Schedule	Assay Dependant
Expected TAT	Assay Dependant
Notes	See individual panel components for more information for those tests.
CPT Code(s)	83519x3, 86255
Lab Section	Reference Lab