## **Acid Phosphate Total**

Order Name: ACID PHS T
Test Number: 2000280
Revision Date: 12/12/2022

TEST NAME			METHODOLOGY	LOINC CODE	
Acid Phosphate To	otal		Enzymatic	12173-1	
SPECIMEN REQU	JIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	1.5 mL (0.8 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen	
Instructions	Allow specimen to clot completely. Transfer 1.5mL of serum into an aliquot tube and freeze.				
	CRITICAL FROZEN. Separate	e specimens must be submitt	ted when multiple tests are ordered.		
	Unacceptable Conditions: Plas	sma. Non-frozen specimens. He	emolyzed specimens.		
	Stability After separation from cells: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 1 month				

GENERAL INFORMATION	
Expected TAT	5 - 8 days
Clinical Use	This test is used in the evaluation of prostate cancer, Padget's disease, hyperthyroidism, and Gaucher's disease.
Performing Labcorp Test Code	829642
Notes	Labcorp Test Code: 829642
CPT Code(s)	84060
Lab Section	Reference Lab

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025