Coffee IgG

Order Name: **COFFEE IGG** Test Number: 5500589 Revision Date: 12/14/2022

TEST NAME			METHODOLOGY	LOINC CODE	
Coffee IgG					
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	0.3 mL (0.1 mL)	Serum	Clot Activator SST	Room Temperature	
Instructions	Specimen Type: Red-top tube or gel-barrier tube				
	Specimen Storage: Room temperature				
	Specimen Collection: Not Av	ailable			
	Special Instructions: Values obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only				
	one assay method be used consistently to monitor each patient's course of therapy.				
	Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available				
GENERAL INFOR	MATION				
Expected TAT	2 - 4 days				

Expected TAT	2 - 4 days
Performing Labcorp Test Code	605134
Notes	Labcorp Test Code: 605134
CPT Code(s)	86001
Lab Section	Reference Lab