

Penicillin Allergy Panel

Order Name: **PENICILIN**
Test Number: 5584380
Revision Date: 08/01/2023

| TEST NAME | METHODOLOGY | LOINC CODE |
|--------------------------|-------------|------------|
| Penicillin Allergy Panel | ImmunoCAP | |

SPECIMEN REQUIREMENTS

| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
|--------------|---|---------------|--------------------|-----------------------|
| Preferred | 1.6 mL (0.8 mL) | Serum | Clot Activator SST | Room Temperature |
| Instructions | Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3) | | | |

GENERAL INFORMATION

| | |
|------------------------------|---|
| Expected TAT | 3-5 Days |
| Performing Labcorp Test Code | 805977, 805999, 602646, 067116 |
| Notes | Testing Performed at Labcorp Phoenix as Individual Allergens. |
| CPT Code(s) | 86003 x 4 |
| Lab Section | Reference Lab |