

EBV (Epstein Barr Virus) Viral Capsid Antibodies IgG and IgM

Order Name: **EBV VCA Ab G/M**

Test Number: 5580926

Revision Date: 04/11/2023

| TEST NAME                  | METHODOLOGY             | LOINC CODE |
|----------------------------|-------------------------|------------|
| EBV VCA IgG                | Chemiluminescence Assay | 5157-3     |
| EBV VCA IgG Interpretation |                         | 24114-1    |
| EBV VCA IgM                | Chemiluminescence Assay | 5159-9     |
| EBV VCA IgM Interpretation |                         | 24115-8    |

| SPECIMEN REQUIREMENTS |  |               |                    |                        |
|-----------------------|--|---------------|--------------------|------------------------|
| Specimen              | Specimen Volume (min)  | Specimen Type | Specimen Container | Transport Environment  |
| Preferred             | 2 mL (0.25mL)  | Serum         | Clot Activator SST | Refrigerated or Frozen |
| Instructions          | Specimen Stability: Room Temperature: Unacceptable, Refrigerated: 7 days, Frozen 30 Days |               |                    |                        |

| GENERAL INFORMATION |   |
|---------------------|---|
| Testing Schedule    | Mon-Sat   |
| Expected TAT        | 1-2 Days  |
| Notes               | EB Virus IgM IFA confirmation testing will no longer be performed due to supply chain availability. The results for this test are not, by themselves, diagnostic and should be considered in association with clinical presentation and clinical history. |
| CPT Code(s)         | 86665x2   |
| Lab Section         | Immunology - Serology   |