

PIVKA-II (Des-gamme-carboxy Prothrombin)

Order Name: **PIVKA-II DCP**

Test Number: 6905109

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
PIVKA-II (Des-gamme-carboxy Prothrombin)	Liquid-phase binding assay (LiBaSys)	34444-0

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.25 mL)	Serum	Clot Activator SST	Frozen
Alternate 1	1 mL (0.25 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	<p>Notes: 0.25 mL (Note: This volume Does NOT allow for repeat testing) Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: FREEZE Specimen Collection: Separate serum from cells and transfer to a plastic transport tube. Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p>			

GENERAL INFORMATION

Expected TAT	5 - 7 days
Performing Labcorp Test Code	141325
Notes	Labcorp Test Code: 141325
CPT Code(s)	83951
Lab Section	Reference Lab