## Labcorp Oklahoma, Inc. Test Directory

## **Streptomycin Antimicrobial Serum Level**

Order Name: Strep Antimicro Srm

Test Number: 6905151
Revision Date: 01/02/2018

TEST NAME METHODOLOGY LOINC CODE

Streptomycin Antimicrobial Serum Level

Quantitative High Performance Liquid Chromatography

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	Specimen Type: Red-top tube Specimen Storage: FREEZE Specimen Collection: Serum should be separated from cells within two hours of venipuncture. Send serum in a plastic transport tube before freezing. Special Instructions: Testing referred to National Jewish Health NJHCO#PKSM Was SMH NOTE: Drug dose and frequency, including time of last dose should be indicated on request form.			

GENERAL INFORMATION		
14-17 days		
815139		
Labcorp Test Code: 815139		
80299		
Reference Lab		

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025