

Neisseria gonorrhoeae PAP - TMA

Order Name: **GC PAP**

Test Number: 6910113

Revision Date: 10/10/2018

TEST NAME	METHODOLOGY	LOINC CODE
Neisseria gonorrhoeae PAP	Transcription-Mediated Amplification	50388-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	PAP specimen	ThinPrep PreservCyt solution	Room Temperature
Instructions	Cervical brush/spatula in ThinPrep Pap test collection kit. ThinPrep PreservCyt solution. <b>Stability:</b> Specimen must be transferred from ThinPrep PAP vial into an Aptima Specimen Transfer tube within 30 days. Once transferred, the specimen is stable in the Aptima vial for 30 days refrigerated.  <b>Please contact your account rep if training and ThinPrep supplies have not been supplied to your practice</b>			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	1-2 days after set-up
CPT Code(s)	87591
Lab Section	Molecular Diagnostic